

## Product datasheet for **TL514329**

### Sept7 Mouse shRNA Plasmid (Locus ID 235072)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Sept7 Mouse shRNA Plasmid (Locus ID 235072)
Locus ID:	235072
Synonyms:	Cdc10; E430034N22
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Sept7 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 235072). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC058587</a> , <a href="#">NM_001205367</a> , <a href="#">NM_009859</a> , <a href="#">NM_001359736</a> , <a href="#">NM_001359737</a> , <a href="#">NM_001359738</a> , <a href="#">NM_009859.1</a> , <a href="#">NM_009859.2</a> , <a href="#">NM_009859.3</a> , <a href="#">NM_009859.4</a> , <a href="#">NM_001205367.1</a>
UniProt ID:	<a href="#">O55131</a>
Summary:	Filament-forming cytoskeletal GTPase. Required for normal organization of the actin cytoskeleton. Required for normal progress through mitosis. Involved in cytokinesis. Required for normal association of CENPE with the kinetochore. Plays a role in ciliogenesis and collective cell movements. Forms a filamentous structure with SEPTIN12, SEPTIN6, SEPTIN2 and probably SEPTIN4 at the sperm annulus which is required for the structural integrity and motility of the sperm tail during postmeiotic differentiation (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).