

## Product datasheet for **TL514665**

### Wnk2 Mouse shRNA Plasmid (Locus ID 75607)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Wnk2 Mouse shRNA Plasmid (Locus ID 75607)
Locus ID:	75607
Synonyms:	1810073P09Rik; AW122246; ESTM15; mKIAA1760; X83337
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Wnk2 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 75607). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001290311</a> , <a href="#">NM_001290313</a> , <a href="#">NM_029361</a> , <a href="#">NM_029361.1</a> , <a href="#">NM_029361.2</a> , <a href="#">NM_029361.3</a> , <a href="#">NM_029361.4</a> , <a href="#">NM_001290313.1</a> , <a href="#">NM_001290311.1</a> , <a href="#">BC046464</a> , <a href="#">BC055795</a> , <a href="#">BC060187</a> , <a href="#">BM941970</a>
UniProt ID:	<a href="#">Q3UH66</a>
Summary:	Serine/threonine kinase which plays an important role in the regulation of electrolyte homeostasis, cell signaling, survival, and proliferation. Acts as an activator and inhibitor of sodium-coupled chloride cotransporters and potassium-coupled chloride cotransporters respectively. Activates SLC12A2, SCNN1A, SCNN1B, SCNN1D and SGK1 and inhibits SLC12A5. Negatively regulates the EGF-induced activation of the ERK/MAPK-pathway and the downstream cell cycle progression. Affects MAPK3/MAPK1 activity by modulating the activity of MAP2K1 and this modulation depends on phosphorylation of MAP2K1 by PAK1. WNK2 acts by interfering with the activity of PAK1 by controlling the balance of the activity of upstream regulators of PAK1 activity, RHOA and RAC1, which display reciprocal activity. [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).