

Product datasheet for **TL514771**

Cd300lf Mouse shRNA Plasmid (Locus ID 246746)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Cd300lf Mouse shRNA Plasmid (Locus ID 246746)
Locus ID:	246746
Synonyms:	CLIM1; CLM-1; CLM1; Digr2; F730004D16Rik; IgSF13; IREM1; LMIR3; Pigr3
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Cd300lf - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 246746). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC057864 , NM_001169153 , NM_001347648 , NM_145634 , NM_145634.1 , NM_145634.2 , NM_145634.3 , NM_001169153.1 , BM227747
UniProt ID:	Q6SIQ7
Summary:	Acts as an inhibitory receptor for myeloid cells and mast cells (PubMed:17438331). Positively regulates the phagocytosis of apoptotic cells (efferocytosis) via phosphatidylserine (PS) recognition; recognizes and binds PS as a ligand which is expressed on the surface of apoptotic cells (PubMed:21865548). Plays an important role in the maintenance of immune homeostasis, by promoting macrophage-mediated efferocytosis and by inhibiting dendritic cell-mediated efferocytosis (PubMed:26768664). Negatively regulates Fc epsilon receptor-dependent mast cell activation and allergic responses via binding to ceramide which acts as a ligand (PubMed:23123064). May act as a coreceptor for interleukin 4 (IL-4). Associates with and regulates IL-4 receptor alpha-mediated responses by augmenting IL-4- and IL-13-induced signaling (PubMed:26124135). Negatively regulates the Toll-like receptor (TLR) signaling mediated by MYD88 and TRIF through activation of PTPN6/SHP-1 and PTPN11/SHP-2 (By similarity). Inhibits osteoclast formation (PubMed:14662855). Induces macrophage cell death upon engagement (PubMed:18097021).[UniProtKB/Swiss-Prot Function]



[View online »](#)

shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

Performance Guaranteed: OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).