

Product datasheet for **TL514855V**

Uchl1 Mouse shRNA Lentiviral Particle (Locus ID 22223)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Uchl1 Mouse shRNA Lentiviral Particle (Locus ID 22223)
Locus ID:	22223
Synonyms:	AW822034; C88048; gad; PGP 9.5; PGP9.5; R75593; UCH-L1; UCHL-1
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Uchl1 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	BC039177 , NM_011670 , NM_011670.1 , NM_011670.2
UniProt ID:	Q9R0P9
Summary:	Ubiquitin-protein hydrolase involved both in the processing of ubiquitin precursors and of ubiquitinated proteins (Probable). This enzyme is a thiol protease that recognizes and hydrolyzes a peptide bond at the C-terminal glycine of ubiquitin (PubMed:12913066). Also binds to free monoubiquitin and may prevent its degradation in lysosomes (PubMed:12913066). The homodimer may have ATP-independent ubiquitin ligase activity (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).