

## Product datasheet for **TL514879**

### Pan2 Mouse shRNA Plasmid (Locus ID 103135)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Pan2 Mouse shRNA Plasmid (Locus ID 103135)
Locus ID:	103135
Synonyms:	1200014O24Rik; AI047843; AW742773; mKIAA0710; Usp52
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Pan2 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 103135). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC075686</a> , <a href="#">BC079841</a> , <a href="#">NM_001252326</a> , <a href="#">NM_001252327</a> , <a href="#">NM_133992</a> , <a href="#">NM_001358822</a> , <a href="#">NM_133992.1</a> , <a href="#">NM_133992.2</a> , <a href="#">NM_133992.3</a> , <a href="#">NM_001252327.1</a> , <a href="#">NM_001252326.1</a> , <a href="#">BC079841.1</a> , <a href="#">BC006064</a> , <a href="#">BC008521</a> , <a href="#">BC020158</a> , <a href="#">BC043659</a>
UniProt ID:	<a href="#">Q8BGF7</a>
Summary:	Catalytic subunit of the poly(A)-nuclease (PAN) deadenylation complex, one of two cytoplasmic mRNA deadenylases involved in general and miRNA-mediated mRNA turnover. PAN specifically shortens poly(A) tails of RNA and the activity is stimulated by poly(A)-binding protein (PABP). PAN deadenylation is followed by rapid degradation of the shortened mRNA tails by the CCR4-NOT complex. Deadenylated mRNAs are then degraded by two alternative mechanisms, namely exosome-mediated 3'-5' exonucleolytic degradation, or deadenylation-dependent mRNA decapping and subsequent 5'-3' exonucleolytic degradation by XRN1 (PubMed:16284618). Also acts as an important regulator of the HIF1A-mediated hypoxic response. Required for HIF1A mRNA stability independent of poly(A) tail length regulation (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).