

Product datasheet for **TL514933**

Gfra4 Mouse shRNA Plasmid (Locus ID 14588)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Gfra4 Mouse shRNA Plasmid (Locus ID 14588)
Locus ID:	14588
Synonyms:	AW060733; G630015H18Rik
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Gfra4 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 14588). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC141874 , NM_001136063 , NM_001271001 , NM_001271002 , NM_020014 , NM_020014.1 , NM_020014.2 , NM_001136063.1 , NM_001136063.2 , NM_001271001.1 , NM_001271002.1 , BC111808 , BC141873
UniProt ID:	Q9JIT2
Summary:	This gene encodes a transmembrane protein that functions as the receptor for persephin, a member of the glial cell line derived neurotrophic factors. The encoded protein undergoes proteolytic processing to generate a glycosylphosphatidylinositol-anchored cell surface coreceptor that forms a complex with the receptor tyrosine kinase Ret. A complete lack of the encoded protein impairs production of thyroid calcitonin and increases the rate of bone formation in young mice. [provided by RefSeq, Jul 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).