

Product datasheet for **TL515021**

Dnaja1 Mouse shRNA Plasmid (Locus ID 15502)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Dnaja1 Mouse shRNA Plasmid (Locus ID 15502)
Locus ID:	15502
Synonyms:	Hsj; HSJ-2; Hsj2; Nedd; Nedd7
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Dnaja1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 15502). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC057876 , NM_001164671 , NM_001164672 , NM_008298 , NM_008298.1 , NM_008298.2 , NM_008298.3 , NM_008298.4 , NM_008298.5 , NM_008298.6 , NM_001164672.1 , NM_001164672.2 , NM_001164671.1 , NM_001164671.2 , BC060653 , BC158024
UniProt ID:	P63037
Summary:	The protein encoded by this gene is a member of the Dnaj family, whose members act as cochaperones of heat shock protein 70. Heat shock proteins facilitate protein folding, trafficking, prevention of aggregation, and proteolytic degradation. Members of this family are characterized by a highly conserved N-terminal J domain, a glycine/phenylalanine-rich region, four CxxCxGxG zinc finger repeats, and a C-terminal substrate-binding domain. The J domain mediates the interaction with heat shock protein 70 to recruit substrates and regulate ATP hydrolysis activity. Mice deficient for this gene display reduced levels of activation-induced deaminase, an enzyme that deaminates deoxycytidine at the immunoglobulin genes during immune responses. In addition, mice lacking this gene exhibit severe defects in spermatogenesis. Several pseudogenes of this gene are found on other chromosomes. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2015]



[View online »](#)

- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).