

## Product datasheet for **TL515070V**

### Hormad1 Mouse shRNA Lentiviral Particle (Locus ID 67981)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Hormad1 Mouse shRNA Lentiviral Particle (Locus ID 67981)
Locus ID:	67981
Synonyms:	4921522K05Rik; Nohma
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Hormad1 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">BC051129</a> , <a href="#">NM_001289532</a> , <a href="#">NM_001289534</a> , <a href="#">NM_001289537</a> , <a href="#">NM_026489</a> , <a href="#">NM_026489.1</a> , <a href="#">NM_026489.2</a> , <a href="#">NM_026489.3</a> , <a href="#">NM_001289534.1</a> , <a href="#">NM_001289537.1</a> , <a href="#">NM_001289532.1</a>
UniProt ID:	<a href="#">Q9D5T7</a>
Summary:	Plays a key role in meiotic progression (PubMed:19686734, PubMed:21079677, PubMed:21478856). Regulates 3 different functions during meiosis: ensures that sufficient numbers of processed DNA double-strand breaks (DSBs) are available for successful homology search by increasing the steady-state numbers of single-stranded DSB ends (PubMed:19686734, PubMed:21079677). Promotes synaptonemal-complex formation independently of its role in homology search (PubMed:19686734, PubMed:21079677). Plays a key role in the male mid-pachytene checkpoint and the female meiotic prophase checkpoint: required for efficient build-up of ATR activity on unsynapsed chromosome regions, a process believed to form the basis of meiotic silencing of unsynapsed chromatin (MSUC) and meiotic prophase quality control in both sexes (PubMed:21478856).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).