

Product datasheet for **TL515129**

Spsb2 Mouse shRNA Plasmid (Locus ID 14794)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Spsb2 Mouse shRNA Plasmid (Locus ID 14794)
Locus ID:	14794
Synonyms:	AI461677; C9; Grcc; Grcc9; SS; SSB2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Spsb2 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 14794). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC002005 , BC010305 , NM_013539 , NM_013539.1 , NM_013539.2
UniProt ID:	O88838
Summary:	This gene encodes a member of the SSB family of proteins that contain a central SPRY (repeats in splA and ryanodine receptors) domain and a C-terminal SOCS (suppressor of cytokine signaling) box. The encoded protein is an adaptor protein in the E3 ubiquitin ligase complex that ubiquitinates inducible nitric oxide synthase and targets it for proteasomal degradation. Mice lacking the encoded protein exhibit lower blood urea nitrogen levels and mild thrombocytopenia due to reduced platelet production. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Apr 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).