

Product datasheet for **TL515143**

Ptrf Mouse shRNA Plasmid (Locus ID 19285)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Ptrf Mouse shRNA Plasmid (Locus ID 19285)
Locus ID:	19285
Synonyms:	2310075E07Rik; AW546441; Cav-p60; Cavin
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Ptrf - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 19285). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC110698 , BC116701 , BC117527 , NM_008986 , NM_001359105 , NM_008986.1 , NM_008986.2 , BC065792 , BC099493
UniProt ID:	O54724
Summary:	Plays an important role in caveolae formation and organization. Essential for the formation of caveolae in all tissues (PubMed:18191225, PubMed:18840361, PubMed:18056712, PubMed:30188967). Core component of the CAVIN complex which is essential for recruitment of the complex to the caveolae in presence of calveolin-1 (CAV1) (PubMed:19546242). Essential for normal oligomerization of CAV1 (PubMed:23652019). Promotes ribosomal transcriptional activity in response to metabolic challenges in the adipocytes and plays an important role in the formation of the ribosomal transcriptional loop (PubMed:27528195). Dissociates transcription complexes paused by DNA-bound TTF1, thereby releasing both RNA polymerase I and pre-RNA from the template (PubMed:9582279, PubMed:11139612). The caveolae biogenesis pathway is required for the secretion of proteins such as GASK1A (PubMed:30188967).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).