

## Product datasheet for **TL515155**

### Dst Mouse shRNA Plasmid (Locus ID 13518)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Dst Mouse shRNA Plasmid (Locus ID 13518)
Locus ID:	13518
Synonyms:	2310001O04Rik; A830042E19Rik; ah; athetoid; AW554249; BP230; Bpag; Bpag1; BPAG1-n; dt; Macf2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Dst - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 13518). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001276764</a> , <a href="#">NM_010081</a> , <a href="#">NM_133833</a> , <a href="#">NM_134448</a> , <a href="#">NM_010081.1</a> , <a href="#">NM_010081.2</a> , <a href="#">BC019791</a> , <a href="#">BC033468</a> , <a href="#">BC048929</a> , <a href="#">BC131678</a> , <a href="#">NM_010081.3</a>
UniProt ID:	<a href="#">Q91ZU6</a>
Summary:	Cytoskeletal linker protein. Acts as an integrator of intermediate filaments, actin and microtubule cytoskeleton networks. Required for anchoring either intermediate filaments to the actin cytoskeleton in neural and muscle cells or keratin-containing intermediate filaments to hemidesmosomes in epithelial cells. The proteins may self-aggregate to form filaments or a two-dimensional mesh. Regulates the organization and stability of the microtubule network of sensory neurons to allow axonal transport. Mediates docking of the dynein/dynactin motor complex to vesicle cargos for retrograde axonal transport through its interaction with TMEM108 and DCTN1.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).