

Product datasheet for **TL515251V**

Smdt1 Mouse shRNA Lentiviral Particle (Locus ID 69029)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Smdt1 Mouse shRNA Lentiviral Particle (Locus ID 69029)
Locus ID:	69029
Synonyms:	1500032L24Rik; Emre
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	1500032L24Rik - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	BC028457 , NM_026914 , NM_026914.1 , BC092006
UniProt ID:	Q9DB10
Summary:	Essential regulatory subunit of the mitochondrial calcium uniporter complex (uniplx), a complex that mediates calcium uptake into mitochondria (PubMed:27001609). Required to bridge the calcium-sensing proteins MICU1 and MICU2 with the calcium-conducting subunit MCU. Plays a central role in regulating the uniplx complex response to intracellular calcium signaling. Acts by mediating activation of MCU and retention of MICU1 to the MCU pore, in order to ensure tight regulation of the uniplx complex and appropriate responses to intracellular calcium signaling (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).