

Product datasheet for **TL515321**

Gramd4 Mouse shRNA Plasmid (Locus ID 223752)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Gramd4 Mouse shRNA Plasmid (Locus ID 223752)
Locus ID:	223752
Synonyms:	9930016O13
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Gramd4 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 223752). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001205353 , NM_172611 , NM_172611.1 , NM_172611.2 , NM_172611.3 , NM_172611.4 , NM_001205353.1 , BC139228 , BC021523 , BC085148 , BC139229
UniProt ID:	Q8CB44
Summary:	Plays a role as a mediator of E2F1-induced apoptosis in the absence of p53/TP53 (By similarity). Inhibits TLR9 response to nucleic acids and regulates TLR9-mediated innate immune response (PubMed:25917084).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .

[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).