

Product datasheet for **TL515445V**

Notch3 Mouse shRNA Lentiviral Particle (Locus ID 18131)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Notch3 Mouse shRNA Lentiviral Particle (Locus ID 18131)
Locus ID:	18131
Synonyms:	AW229011; hpbk; N3
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Notch3 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_008716 , NM_008716.1 , NM_008716.2 , BC172738 , NM_008716.3
UniProt ID:	Q61982
Summary:	Functions as a receptor for membrane-bound ligands Jagged1, Jagged2 and Delta1 to regulate cell-fate determination. Upon ligand activation through the released notch intracellular domain (NICD) it forms a transcriptional activator complex with RBPJ/RBPSUH and activates genes of the enhancer of split locus. Affects the implementation of differentiation, proliferation and apoptotic programs (By similarity). May play a role during CNS development.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:

GFP signal was observed under microscope at 48 hours after transduction of TL515445B virus into HEK293 cells. TL515445B virus was prepared using lenti-shRNA TL515445B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL515445D] virus into HEK293 cells. [TL515445D] virus was prepared using lenti-shRNA [TL515445D] and [TR30037] packaging kit.