

Product datasheet for **TL515601**

Cldn10 Mouse shRNA Plasmid (Locus ID 58187)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Cldn10 Mouse shRNA Plasmid (Locus ID 58187)
Locus ID:	58187
Synonyms:	6720456116Rik; Cldn; Cldn1; Cldn10a; Cldn10b; D14Ertd728; D14Ertd728e
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Cldn10 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 58187). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC021770 , BC029019 , NM_001160096 , NM_001160097 , NM_001160098 , NM_001160099 , NM_021386 , NM_023878 , NM_001160097.1 , NM_021386.1 , NM_021386.2 , NM_021386.3 , NM_021386.4 , NM_023878.1 , NM_023878.2 , NM_001160098.1 , NM_001160099.1 , BC012238
UniProt ID:	Q9Z0S6
Summary:	This intronless gene encodes a member of the claudin family. Claudins are integral membrane proteins and components of tight junction strands. Tight junction strands serve as a physical barrier to prevent solutes and water from passing freely through the paracellular space between epithelial or endothelial cell sheets, and also play critical roles in maintaining cell polarity and signal transductions. Six alternatively spliced transcript variants have been identified, which encode different isoforms with distinct electric charge of the first extracellular loop and with or without the fourth transmembrane region. These isoforms exhibit distinct localization and function in paracellular anion or cation permeability. [provided by RefSeq, Aug 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).