

Product datasheet for **TL515609V**

Mdm4 Mouse shRNA Lentiviral Particle (Locus ID 17248)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Mdm4 Mouse shRNA Lentiviral Particle (Locus ID 17248)
Locus ID:	17248
Synonyms:	4933417N07Rik; AA414968; AL023055; AU018793; AU021806; C85810; Mdmx
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Mdm4 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	BC003750 , NM_001302801 , NM_001302802 , NM_001302803 , NM_001302804 , NM_008575 , NR_126506 , NM_008575.1 , NM_008575.2 , NM_008575.3 , NM_008575.4 , NM_001302804.1 , NM_001302803.1 , NM_001302801.1 , NM_001302802.1 , BC023912 , BC025868 , BC058680
UniProt ID:	O35618
Summary:	This gene encodes a protein that has been shown to negatively regulate the activity of the tumor suppressor protein p53. Homozygous knockout mice exhibit embryonic lethality as a result of p53-dependent apoptosis and cell cycle arrest. Amplification of this gene or overexpression of the encoded protein has been linked to a range of human cancers. A pseudogene has been identified on the X chromosome. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Nov 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).