

## Product datasheet for **TL515631**

### **Fbn1 Mouse shRNA Plasmid (Locus ID 14118)**

#### **Product data:**

<b>Product Type:</b>	shRNA Plasmids
<b>Product Name:</b>	Fbn1 Mouse shRNA Plasmid (Locus ID 14118)
<b>Locus ID:</b>	14118
<b>Synonyms:</b>	AI536462; B430209H23; Fib-; Fib-1; Tsk
<b>Vector:</b>	pGFP-C-shLenti (TR30023)
<b>E. coli Selection:</b>	Chloramphenicol (34 ug/ml)
<b>Mammalian Cell Selection:</b>	Puromycin
<b>Format:</b>	Lentiviral plasmids
<b>Components:</b>	Fbn1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 14118). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
<b>RefSeq:</b>	<a href="#">NM_007993</a> , <a href="#">NM_007993.1</a> , <a href="#">NM_007993.2</a> , <a href="#">BC003905</a> , <a href="#">BC055932</a>
<b>UniProt ID:</b>	<a href="#">Q61554</a>
<b>Summary:</b>	This gene encodes a member of the fibrillin family of proteins. The encoded preproprotein is proteolytically processed to generate two proteins including the extracellular matrix component fibrillin-1 and the protein hormone asprosin. Fibrillin-1 is an extracellular matrix glycoprotein that serves as a structural component of calcium-binding microfibrils. Asprosin, secreted by white adipose tissue, has been shown to regulate glucose homeostasis. Homozygous knockout mice for this gene exhibit impaired aortic development and early postnatal death, which was attributed to a deficiency in the fibrillin-1 protein. Mice with a hypomorphic allele of this gene exhibit impaired glucose homeostasis, likely due to a reduction in serum asprosin levels. [provided by RefSeq, Apr 2016]
<b>shRNA Design:</b>	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .




[View online »](#)


**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

**Product images:**


GFP signal was observed under microscope at 48 hours after transduction of TL515631A virus into HEK293 cells. TL515631A virus was prepared using lenti-shRNA TL515631A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL515631B virus into HEK293 cells. TL515631B virus was prepared using lenti-shRNA TL515631B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL515631C] virus into HEK293 cells. [TL515631C] virus was prepared using lenti-shRNA [TL515631C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL515631D] virus into HEK293 cells. [TL515631D] virus was prepared using lenti-shRNA [TL515631D] and [TR30037] packaging kit.