

## Product datasheet for **TL515656**

### Mag Mouse shRNA Plasmid (Locus ID 17136)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Mag Mouse shRNA Plasmid (Locus ID 17136)
Locus ID:	17136
Synonyms:	Gm; Gma; sigle; siglec-4a
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Mag - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 17136). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC101944</a> , <a href="#">BC103521</a> , <a href="#">NM_010758</a> , <a href="#">NM_010758.1</a> , <a href="#">NM_010758.2</a> , <a href="#">NM_010758.3</a> , <a href="#">BC103522</a> , <a href="#">BC103523</a>
UniProt ID:	<a href="#">P20917</a>
Summary:	This gene encodes a type I membrane protein and member of the immunoglobulin-like superfamily. It is expressed in myelinating glial cells, including oligodendrocytes of the central nervous system and Schwann cells of the peripheral nervous system. Mice lacking the encoded protein express abundant myelin, but suffer long-term axon degeneration, altered distribution of channels and adhesion molecules at nodes of Ranvier, and altered axon cytoskeletal structure. While not required for myelination, the encoded protein enhances axon-myelin stability, helps to structure nodes of Ranvier, and regulates the axon cytoskeleton. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).