

## Product datasheet for **TL515657**

### Atxn3 Mouse shRNA Plasmid (Locus ID 110616)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Atxn3 Mouse shRNA Plasmid (Locus ID 110616)
Locus ID:	110616
Synonyms:	2210008M02Rik; AI463012; AI647473; ataxin-3; ATX3; Mjd; MJD1; Sca3
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Atxn3 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 110616). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC087880</a> , <a href="#">NM_001167914</a> , <a href="#">NM_029705</a> , <a href="#">NM_001167914.1</a> , <a href="#">NM_029705.1</a> , <a href="#">NM_029705.2</a> , <a href="#">NM_029705.3</a> , <a href="#">BC049743</a>
UniProt ID:	<a href="#">Q9CVD2</a>
Summary:	Deubiquitinating enzyme involved in protein homeostasis maintenance, transcription, cytoskeleton regulation, myogenesis and degradation of misfolded chaperone substrates (By similarity). Binds long polyubiquitin chains and trims them, while it has weak or no activity against chains of 4 or less ubiquitins (By similarity). Involved in degradation of misfolded chaperone substrates via its interaction with STUB1/CHIP: recruited to monoubiquitinated STUB1/CHIP, and restricts the length of ubiquitin chain attached to STUB1/CHIP substrates and preventing further chain extension (PubMed:21855799). Interacts with key regulators of transcription and represses transcription: acts as a histone-binding protein that regulates transcription (By similarity). Regulates autophagy via the deubiquitination of 'Lys-402' of BECN1 leading to the stabilization of BECN1 (PubMed:28445460).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).