

Product datasheet for **TL515851**

Dmd Mouse shRNA Plasmid (Locus ID 13405)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Dmd Mouse shRNA Plasmid (Locus ID 13405)
Locus ID:	13405
Synonyms:	Dp7; Dp42; Dp71; Dp427; DXSmh7; DXSmh9; dy; dys; mdx; pk; pke
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Dmd - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 13405). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC085236
Summary:	This gene encodes a large, rod-like cytoskeletal protein which is found at the inner surface of muscle fibers in skeletal and cardiac muscles. The encoded protein, dystrophin, is part of the dystrophin-glycoprotein complex, which bridges the inner cytoskeleton (F-actin) and the extra-cellular matrix. This protein is required for proper development and organization of myofibers as contractile units in striated muscles. Mutations in the human gene cause Duchenne and Becker Muscular Dystrophies and a form of heart disease called DMD-associated dilated cardiomyopathy. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined. [provided by RefSeq, Sep 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).