

Product datasheet for **TL515918**

Cryga Mouse shRNA Plasmid (Locus ID 12964)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Cryga Mouse shRNA Plasmid (Locus ID 12964)
Locus ID:	12964
Synonyms:	Cryg-; Cryg-4; DGcry-; DGcry-4; Sec; Secc
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Cryga - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 12964). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC056430 , NM_007774 , NM_007774.1 , NM_007774.2 , NM_007774.3
UniProt ID:	P04345
Summary:	Three main families of major soluble proteins, the alpha, beta and gamma crystallins, are ubiquitously expressed in vertebrate lenses. This gene encodes a member of the gamma-crystallin family of proteins which may function as a structural component of the eye lens. Gamma-crystallins are a homogeneous group of highly symmetrical, monomeric proteins typically lacking connecting peptides and terminal extensions. They are differentially regulated after early development. Five gamma-crystallin genes (gamma-A through gamma-E) are tandemly organized in a genomic segment as a gene cluster in the mouse. Another gamma-crystallin gene (gamma-F) is found some distance upstream of the cluster on the same chromosome. Whether due to aging or mutations in specific genes, gamma-crystallins have been involved in cataract formation. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).