

Product datasheet for **TL515957V**

Sult2b1 Mouse shRNA Lentiviral Particle (Locus ID 54200)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Sult2b1 Mouse shRNA Lentiviral Particle (Locus ID 54200)
Locus ID:	54200
Synonyms:	AI326997; BB173635; ST2B1; SULT2B
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Sult2b1 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	BC009811 , BC009813 , NM_017465 , NM_001360784 , NM_017465.1 , NM_017465.2 , BM899747 , NM_017465.3
UniProt ID:	O35400
Summary:	Sulfotransferase that utilizes 3'-phospho-5'-adenylyl sulfate (PAPS) as sulfonate donor to catalyze the sulfate conjugation of many hormones, neurotransmitters, drugs and xenobiotic compounds. Sulfonation increases the water solubility of most compounds, and therefore their renal excretion, but it can also result in bioactivation to form active metabolites. Sulfates hydroxysteroids such as dehydroepiandrosterone. Isoform 1 is required for production of cholesterol sulfate essential for normal skin development whereas isoform 2 produces pregnenolone sulfate, an essential neurosteroid during development of the central nervous system. Plays a role in epidermal cholesterol metabolism and in the regulation of epidermal proliferation and differentiation.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).