

## Product datasheet for **TL515966**

### Nhej1 Mouse shRNA Plasmid (Locus ID 75570)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Nhej1 Mouse shRNA Plasmid (Locus ID 75570)
Locus ID:	75570
Synonyms:	1700029B21Rik; cernunnos; XLF
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Nhej1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 75570). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC006063</a> , <a href="#">BC107252</a> , <a href="#">BC107253</a> , <a href="#">BC132359</a> , <a href="#">BC132361</a> , <a href="#">NM_029342</a> , <a href="#">NM_029342.1</a> , <a href="#">NM_029342.2</a> , <a href="#">NM_029342.3</a> , <a href="#">NM_029342.4</a>
UniProt ID:	<a href="#">Q3KNJ2</a>
Summary:	DNA repair protein involved in DNA nonhomologous end joining (NHEJ) required for double-strand break (DSB) repair and V(D)J recombination. May serve as a bridge between XRCC4 and the other NHEJ factors located at DNA ends, or may participate in reconfiguration of the end bound NHEJ factors to allow XRCC4 access to the DNA termini. It may act in concert with XRCC6/XRCC5 (Ku) to stimulate XRCC4-mediated joining of blunt ends and several types of mismatched ends that are noncomplementary or partially complementary (PubMed:17360556). Binds DNA in a length-dependent manner (By similarity). [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).