

Product datasheet for **TL515996**

Abhd5 Mouse shRNA Plasmid (Locus ID 67469)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Abhd5 Mouse shRNA Plasmid (Locus ID 67469)
Locus ID:	67469
Synonyms:	1300003D03Rik; 2010002J10Rik; CDS; CGI-58; IECN5; NCIE2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Abhd5 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 67469). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC037063 , NM_026179 , NM_001359207 , NR_152865 , NM_026179.1 , NM_026179.2 , BC006683 , BC080714
UniProt ID:	Q9DBL9
Summary:	Coenzyme A-dependent lysophosphatidic acid acyltransferase that catalyzes the transfert of an acyl group on a lysophosphatidic acid (PubMed:19801371). Functions preferentially with 1-oleoyl-lysophosphatidic acid followed by 1-palmitoyl-lysophosphatidic acid, 1-stearoyl-lysophosphatidic acid and 1-arachidonoyl-lysophosphatidic acid as lipid acceptor (PubMed:19801371). Functions preferentially with arachidonoyl-CoA followed by oleoyl-CoA as acyl group donors (PubMed:19801371). Functions in phosphatidic acid biosynthesis (By similarity). May regulate the cellular storage of triacylglycerol through activation of the phospholipase PNPLA2 (PubMed:16679289). Involved in keratinocyte differentiation (By similarity). Regulates lipid droplet fusion (PubMed:26083785).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).