

Product datasheet for **TL516111**

Wnt9b Mouse shRNA Plasmid (Locus ID 22412)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Wnt9b Mouse shRNA Plasmid (Locus ID 22412)
Locus ID:	22412
Synonyms:	clf; clf1; wnt-14b; wnt-15; Wnt14b; Wnt15
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Wnt9b - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 22412). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC110482 , NM_011719 , NM_011719.1 , NM_011719.2 , NM_011719.3 , NM_011719.4
UniProt ID:	O35468
Summary:	Ligand for members of the frizzled family of seven transmembrane receptors (Probable). Functions in the canonical Wnt/beta-catenin signaling pathway (PubMed:22461561, PubMed:16054034, PubMed:17537789). Required for normal embryonic kidney development, and for normal development of the urogenital tract, including uterus and part of the oviduct and the upper vagina in females, and epididymis and vas deferens in males (PubMed:16054034). Activates a signaling cascade in the metanephric mesenchyme that induces tubulogenesis (PubMed:16054034, PubMed:17537789). Acts upstream of WNT4 in the signaling pathways that mediate development of kidney tubules and the Mullerian ducts (PubMed:16054034). Plays a role in cranofacial development and is required for normal fusion of the palate during embryonic development (PubMed:16054034, PubMed:22461561, PubMed:25257647).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).