

## Product datasheet for **TL516130**

### Gskip Mouse shRNA Plasmid (Locus ID 66787)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Gskip Mouse shRNA Plasmid (Locus ID 66787)
Locus ID:	66787
Synonyms:	4933433P14Rik
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Gskip - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 66787). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC031494</a> , <a href="#">NM_178613</a> , <a href="#">NM_178613.1</a> , <a href="#">NM_178613.2</a> , <a href="#">NM_178613.3</a>
UniProt ID:	<a href="#">Q8BGR8</a>
Summary:	A-kinase anchoring protein for GSK3B and PKA that regulates or facilitates their kinase activity towards their targets. The ternary complex enhances Wnt-induced signaling by facilitating the GSK3B- and PKA-induced phosphorylation of beta-catenin leading to beta-catenin degradation and stabilization respectively. Upon cAMP activation, the ternary complex contributes to neuroprotection against oxidative stress-induced apoptosis by facilitating the PKA-induced phosphorylation of DML1 and PKA-induced inactivation of GSK3B. During neurite outgrowth promotes neuron proliferation; while increases beta-catenin-induced transcriptional activity through GSK3B kinase activity inhibition, reduces N-cadherin level to promote cell cycle progression (By similarity). May play a role in cleft palate formation and is required for postnatal life through modulation of the activity of GSK3B during development (PubMed:26582204).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).