

Product datasheet for **TL516140**

Cd177 Mouse shRNA Plasmid (Locus ID 68891)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Cd177 Mouse shRNA Plasmid (Locus ID 68891)
Locus ID:	68891
Synonyms:	1190003K14Rik; Pdp3
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Cd177 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 68891). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC027283 , NM_026862 , NM_026862.1 , NM_026862.3 , BC005615
UniProt ID:	Q8R2S8
Summary:	In association with beta-2 integrin heterodimer ITGAM/CD11b and ITGB2/CD18, mediates activation of TNF-alpha primed neutrophils including degranulation and superoxide production (By similarity). In addition, by preventing beta-2 integrin internalization and attenuating chemokine signaling favors adhesion over migration (By similarity). Heterophilic interaction with PECAM1 on endothelial cells plays a role in neutrophil transendothelial migration in vitro (By similarity). However, appears to be dispensable for neutrophil recruitment caused by bacterial infection in vivo (PubMed:25359465). Acts as a receptor for the mature form of protease PRTN3 allowing its display at the cell surface of neutrophils (By similarity). By displaying PRTN3 at the neutrophil cell surface, may play a role in enhancing endothelial cell junctional integrity and thus vascular integrity during neutrophil diapedesis (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).