

## Product datasheet for **TL516592**

### Wwp2 Mouse shRNA Plasmid (Locus ID 66894)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Wwp2 Mouse shRNA Plasmid (Locus ID 66894)
Locus ID:	66894
Synonyms:	1300010O06Rik; AA690238; AIP2; AW554328
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Wwp2 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 66894). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC039921</a> , <a href="#">BC048184</a> , <a href="#">NM_025830</a> , <a href="#">NM_025830.1</a> , <a href="#">NM_025830.2</a> , <a href="#">NM_025830.3</a> , <a href="#">BC004712</a> , <a href="#">BC018496</a> , <a href="#">BC021524</a> , <a href="#">NM_025830.4</a>
UniProt ID:	<a href="#">Q9DBH0</a>
Summary:	E3 ubiquitin-protein ligase which accepts ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates. Polyubiquitinates POU5F1 by 'Lys-63'-linked conjugation and promotes it to proteasomal degradation; regulates POU5F1 protein level during differentiation of embryonal carcinoma cells (ECCs) but not in undifferentiated ECCs and embryonic stem cells (ESCs). Ubiquitinates EGR2 and promotes it to proteasomal degradation; in T-cells the ubiquitination inhibits activation-induced cell death. Ubiquitinates SLC11A2; the ubiquitination is enhanced by presence of NDFIP1 and NDFIP2. Ubiquitinates RPB1 and promotes it to proteasomal degradation.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).