

Product datasheet for **TL516754**

Tjp1 Mouse shRNA Plasmid (Locus ID 21872)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Tjp1 Mouse shRNA Plasmid (Locus ID 21872)
Locus ID:	21872
Synonyms:	ZO1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Tjp1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 21872). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001163574 , NM_009386 , NM_001163574.1 , NM_009386.1 , NM_009386.2 , BC138028 , BC013769 , BC145100
UniProt ID:	P39447
Summary:	Tjp1, TjpP2, and Tjp3 are closely related scaffolding proteins that link tight junction (TJ) transmembrane proteins such as claudins, junctional adhesion molecules, and occludin to the actin cytoskeleton (By similarity). The tight junction acts to limit movement of substances through the paracellular space and as a boundary between the compositionally distinct apical and basolateral plasma membrane domains of epithelial and endothelial cells. Necessary for lumenogenesis, and particularly efficient epithelial polarization and barrier formation (By similarity). Plays a role in the regulation of cell migration by targeting Cdc42bpb to the leading edge of migrating cells (By similarity). Plays an important role in podosome formation and associated function, thus regulating cell adhesion and matrix remodeling (By similarity). With Tjp2 and Tjpp3, participates to the junctional retention and stability of the transcription factor Dbpa, but is not involved in its shuttling to the nucleus (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).