

## Product datasheet for **TL516910**

### Lgals8 Mouse shRNA Plasmid (Locus ID 56048)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Lgals8 Mouse shRNA Plasmid (Locus ID 56048)
Locus ID:	56048
Synonyms:	1200015E08Rik; AI326142; D13Ert524e; Lgals-8
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Lgals8 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 56048). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC040243</a> , <a href="#">NM_001199043</a> , <a href="#">NM_001291055</a> , <a href="#">NM_001291057</a> , <a href="#">NM_001291060</a> , <a href="#">NM_018886</a> , <a href="#">NM_018886.1</a> , <a href="#">NM_018886.2</a> , <a href="#">NM_018886.3</a> , <a href="#">NM_018886.4</a> , <a href="#">NM_018886.5</a> , <a href="#">NM_001291057.1</a> , <a href="#">NM_001291060.1</a> , <a href="#">NM_001199043.1</a> , <a href="#">NM_001291055.1</a>
UniProt ID:	<a href="#">Q9JL15</a>
Summary:	Beta-galactoside-binding lectin that acts as a sensor of membrane damage caused by infection and restricts the proliferation of infecting pathogens by targeting them for autophagy. Detects membrane rupture by binding beta-galactoside ligands located on the luminal side of the endosome membrane; these ligands becoming exposed to the cytoplasm following rupture. Restricts infection by initiating autophagy via interaction with CALCOCO2/NDP52. Required to restrict infection of bacterial invasion such as <i>S.typhimurium</i> . Also required to restrict infection of Picornaviridae viruses. Has a marked preference for 3'-O-sialylated and 3'-O-sulfated glycans.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).