

Product datasheet for **TL516973V**

Bace2 Mouse shRNA Lentiviral Particle (Locus ID 56175)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Bace2 Mouse shRNA Lentiviral Particle (Locus ID 56175)
Locus ID:	56175
Synonyms:	1110059C24Rik; A; AEPLC; AI850424; AL; ALP56; ARP1; AS; ASP1; ASP21; BAE; BAE2; CDA1; CDA13; CEA; CEAP1; DRAP
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Bace2 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	BC120773 , BC131947 , NM_019517 , NM_019517.1 , NM_019517.2 , NM_019517.3 , NM_019517.4 , NM_019517.5
UniProt ID:	Q9JL18
Summary:	This gene encodes a member of the peptidase A1 family of aspartic proteases. The encoded preproprotein undergoes proteolytic processing to generate an active endopeptidase enzyme. This transmembrane protease catalyzes the proteolysis of amyloid precursor protein to produce amyloid beta peptide. Mice lacking the encoded product exhibit increased pancreatic beta cell mass and improved glucose tolerance due to increased insulin secretion. [provided by RefSeq, Jul 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).