

Product datasheet for **TL516978V**

Crh Mouse shRNA Lentiviral Particle (Locus ID 12918)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Crh Mouse shRNA Lentiviral Particle (Locus ID 12918)
Locus ID:	12918
Synonyms:	CR; CRF; Gm1347
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Crh - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	BC116921 , BC119036 , NM_205769 , NM_205769.1 , NM_205769.2 , NM_205769.3
UniProt ID:	Q8CIT0
Summary:	This gene encodes a member of the corticotropin-releasing factor family and preproprotein that is proteolytically processed to generate a mature protein product. This protein product is a neuropeptide hormone that binds to the corticotropin releasing hormone receptors (CRHR1 and CRHR2) to stimulate the release of adrenocorticotrophic hormone from the pituitary gland in response to stress. The encoded protein may also regulate angiogenesis and inflammation. Homozygous knockout mice for this gene exhibit reduced corticosterone levels while the offspring of these mice die perinatally. [provided by RefSeq, Aug 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).