

## Product datasheet for **TL516988**

### Trpm7 Mouse shRNA Plasmid (Locus ID 58800)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Trpm7 Mouse shRNA Plasmid (Locus ID 58800)
Locus ID:	58800
Synonyms:	2310022G15Rik; 4833414K03Rik; 5033407O22Rik; CHAK; CHAK1; Ltp7; LTrpC-7; Ltrpc7
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Trpm7 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 58800). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001164325</a> , <a href="#">NM_021450</a> , <a href="#">NM_001164325.1</a> , <a href="#">NM_021450.1</a> , <a href="#">NM_021450.2</a> , <a href="#">BC150709</a> , <a href="#">BC009137</a> , <a href="#">BC026440</a> , <a href="#">BC145122</a>
UniProt ID:	<a href="#">Q923J1</a>
Summary:	Essential ion channel and serine/threonine-protein kinase. Divalent cation channel permeable to calcium and magnesium. Has a central role in magnesium ion homeostasis and in the regulation of anoxic neuronal cell death. Involved in TNF-induced necroptosis downstream of MLKL by mediating calcium influx. The kinase activity is essential for the channel function. May be involved in a fundamental process that adjusts plasma membrane divalent cation fluxes according to the metabolic state of the cell. Phosphorylates annexin A1 (ANXA1). [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).