

Product datasheet for **TL517228V**

Hnrnpm Mouse shRNA Lentiviral Particle (Locus ID 76936)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Hnrnpm Mouse shRNA Lentiviral Particle (Locus ID 76936)
Locus ID:	76936
Synonyms:	2610023M21Rik; AA409009; Hnrpm; mKIAA4193
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Hnrnpm - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	BC065172 , NM_001109913 , NM_029804 , NM_029804.1 , NM_029804.2 , NM_029804.3 , NM_001109913.1 , BC065172.1 , BC005758 , BC017539 , BC035525 , BC054785
UniProt ID:	Q9D0E1
Summary:	Pre-mRNA binding protein in vivo, binds avidly to poly(G) and poly(U) RNA homopolymers in vitro. Involved in splicing. Acts as a receptor for carcinoembryonic antigen in Kupffer cells, may initiate a series of signaling events leading to tyrosine phosphorylation of proteins and induction of IL-1 alpha, IL-6, IL-10 and tumor necrosis factor alpha cytokines (By similarity). [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).