

## Product datasheet for **TL517332V**

### Cyp2f2 Mouse shRNA Lentiviral Particle (Locus ID 13107)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Cyp2f2 Mouse shRNA Lentiviral Particle (Locus ID 13107)
Locus ID:	13107
Synonyms:	Cyp2f
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Cyp2f2 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">BC011089</a> , <a href="#">BC024742</a> , <a href="#">NM_007817</a> , <a href="#">NM_007817.1</a> , <a href="#">NM_007817.2</a>
UniProt ID:	<a href="#">P33267</a>
Summary:	Involved in the regio- and stereoselective transformation of naphthalene to trans-1R-hydroxy-2R-glutathionyl-1,2-dihydronaphthalene in the presence of glutathione and glutathione S-transferases. It specifically catalyzes the production of a very reactive and potentially toxic intermediate, the 2R,2S arene oxide, that is associated with necrosis of the unciliated bronchiolar epithelial cells or Clara cells in lung.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).