

Product datasheet for TL517565V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Thoc2 Mouse shRNA Lentiviral Particle (Locus ID 331401)

Product data:

Product Type: shRNA Lentiviral Particles

Product Name: Thoc2 Mouse shRNA Lentiviral Particle (Locus ID 331401)

Locus ID: 331401

Synonyms: 6330441O12Rik; D130005M13Rik; Gm1139; Gm1793; Tho2

Vector: pGFP-C-shLenti (TR30023)

Format: Lentiviral particles

Components: Thoc2 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble

control), 0.5 ml each, >10^7 TU/ml.

RefSeq: <u>NM 001033422</u>, <u>NM 001033422.1</u>

UniProt ID: B1AZI6

Summary: Required for efficient export of polyadenylated RNA and spliced mRNA. Acts as component of

the THO subcomplex of the TREX complex which is thought to couple mRNA transcription, processing and nuclear export, and which specifically associates with spliced mRNA and not with unspliced pre-mRNA. TREX is recruited to spliced mRNAs by a transcription-independent mechanism, binds to mRNA upstream of the exon-junction complex (EJC) and is recruited in a

splicing- and cap-dependent manner to a region near the 5' end of the mRNA where it

functions in mRNA export to the cytoplasm via the TAP/NFX1 pathway. Plays a role for proper

neuronal development.[UniProtKB/Swiss-Prot Function]

shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To

be certain that your variant of interest is targeted, please contact <u>techsupport@origene.com</u>. If you need a special design or shRNA sequence, please utilize our <u>custom shRNA service</u>.







Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).