

Product datasheet for **TL517682**

Agbl4 Mouse shRNA Plasmid (Locus ID 78933)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Agbl4 Mouse shRNA Plasmid (Locus ID 78933)
Locus ID:	78933
Synonyms:	4930578N11Rik; 4931433A01Rik
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Agbl4 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 78933). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001048189 , NM_001284190 , NM_030231 , NM_030231.1 , NM_030231.2 , NM_030231.3 , NM_001048189.1 , NM_001048189.2 , NM_001048189.3 , NM_001048189.4 , NM_001284190.1 , NM_001284190.2
UniProt ID:	Q09LZ8
Summary:	Metallocarboxypeptidase that mediates deglutamylation of target proteins (PubMed:17244818, PubMed:21074048, PubMed:25103237, PubMed:26829768). Catalyzes the deglutamylation of polyglutamate side chains generated by post-translational polyglutamylation in proteins such as tubulins (PubMed:17244818). Also removes polyglutamates from the carboxy-terminus of target proteins such as MYLK (PubMed:21074048). Mediates deglutamylation of CGAS, regulating the antiviral activity of CGAS (PubMed:26829768). Acts as a long-chain deglutamylase and specifically shortens long polyglutamate chains, while it is not able to remove the branching point glutamate, a process catalyzed by AGBL5/CCP5 (PubMed:25103237).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).