

Product datasheet for **TL517918**

Pacs2 Mouse shRNA Plasmid (Locus ID 217893)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Pacs2 Mouse shRNA Plasmid (Locus ID 217893)
Locus ID:	217893
Synonyms:	6720425G15Rik; AW208651; AW489976; E230011J22; mKIAA0602; Pacs11
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Pacs2 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 217893). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001081170 , NM_001291444 , NM_001291445 , NM_001081170.1 , NM_001291445.1 , NM_001291444.1 , BC026674 , BC043302 , BC051499 , BC059076 , NM_198165
UniProt ID:	Q3V3Q7
Summary:	Multifunctional sorting protein that controls the endoplasmic reticulum (ER)-mitochondria communication, including the apposition of mitochondria with the ER and ER homeostasis. In addition, in response to apoptotic inducer, translocates BIB to mitochondria, which initiates a sequence of events including the formation of mitochondrial truncated BID, the release of cytochrome c, the activation of caspase-3 thereby causing cell death. May also involved in ion channel trafficking, directing acidic cluster-containing ion channels to distinct subcellular compartments (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).