

Product datasheet for **TL517927**

Heatr5b Mouse shRNA Plasmid (Locus ID 320473)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Heatr5b Mouse shRNA Plasmid (Locus ID 320473)
Locus ID:	320473
Synonyms:	2010013B10Rik; A230048G03Rik; AV071430; D330050P16Rik; mKIAA1414
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Heatr5b - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 320473). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC019508 , BC086652 , NM_001081179 , NM_001081179.1 , NM_028084 , NM_177168 , NM_177332
UniProt ID:	Q8C547
Summary:	Component of clathrin-coated vesicles (By similarity). Component of the aftiphilin/p200/gamma-synergic complex, which plays roles in AP1G1/AP-1-mediated protein trafficking including the trafficking of transferrin from early to recycling endosomes, and the membrane trafficking of furin and the lysosomal enzyme cathepsin D between the trans-Golgi network (TGN) and endosomes (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).