

Product datasheet for **TL517943**

Arid4a Mouse shRNA Plasmid (Locus ID 238247)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Arid4a Mouse shRNA Plasmid (Locus ID 238247)
Locus ID:	238247
Synonyms:	A630009N03; A630067N03Rik; MmRBBP1; Rbbp1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Arid4a - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 238247). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001081195 , NM_001081195.1 , BC157975 , BC025436 , BC037715 , BC058209 , BC080780 , BC158109 , BC160195
UniProt ID:	F8VPQ2
Summary:	DNA-binding protein which modulates activity of several transcription factors including RB1 (retinoblastoma-associated protein) and AR (androgen receptor) (PubMed:17043311, PubMed:23487765). May function as part of an mSin3A repressor complex (By similarity). Has no intrinsic transcriptional activity (PubMed:23487765). Plays a role in the regulation of epigenetic modifications at the PWS/AS imprinting center near the SNRPN promoter, where it might function as part of a complex with RB1 and ARID4B (PubMed:17043311). Involved in spermatogenesis, together with ARID4B, where it acts as a transcriptional coactivator for AR and enhances expression of genes required for sperm maturation (PubMed:23487765). Regulates expression of the tight junction protein CLDN3 in the testis, which is important for integrity of the blood-testis barrier (PubMed:23487765). Plays a role in myeloid homeostasis where it regulates the histone methylation state of bone marrow cells and expression of various genes involved in hematopoiesis (PubMed:18728284). May function as a leukemia suppressor (PubMed:18728284).[UniProtKB/Swiss-Prot Function]



[View online »](#)

shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

Performance Guaranteed: OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).