

Product datasheet for **TL517968**

Erc6 Mouse shRNA Plasmid (Locus ID 319955)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Erc6 Mouse shRNA Plasmid (Locus ID 319955)
Locus ID:	319955
Synonyms:	4732403I04; C130058G22Rik; CSB
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Erc6 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 319955). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001081221 , NM_001081221.1 , BC132447 , BC167234 , NM_177043 , NM_001081221.2
UniProt ID:	F8VPZ5



[View online »](#)

Summary:

Essential factor involved in transcription-coupled nucleotide excision repair which allows RNA polymerase II-blocking lesions to be rapidly removed from the transcribed strand of active genes (By similarity). Upon DNA-binding, it locally modifies DNA conformation by wrapping the DNA around itself, thereby modifying the interface between stalled RNA polymerase II and DNA (By similarity). It is required for transcription-coupled repair complex formation. It recruits the CSA complex (DCX(ERCC8) complex), nucleotide excision repair proteins and EP300 to the sites of RNA polymerase II-blocking lesions (By similarity). Plays an important role in regulating the choice of the DNA double-strand breaks (DSBs) repair pathway and G2/M checkpoint activation; DNA-dependent ATPase activity is essential for this function (By similarity). Regulates the DNA repair pathway choice by inhibiting non-homologous end joining (NHEJ), thereby promoting the homologous recombination (HR)-mediated repair of DSBs during the S/G2 phases of the cell cycle (By similarity). Mediates the activation of the ATM- and CHEK2-dependent DNA damage responses thus preventing the premature exit from the G2/M checkpoint (By similarity). Acts as a chromatin remodeler at DSBs; DNA-dependent ATPase-dependent activity is essential for this function (By similarity). Remodels chromatin by evicting histones from chromatin flanking DSBs, limiting RIF1 accumulation at DSBs thereby promoting BRCA1-mediated HR (By similarity). Required for stable recruitment of ELOA and CUL5 to DNA damage sites (By similarity). Involved in UV-induced translocation of ERCC8 to the nuclear matrix (By similarity). Essential for neuronal differentiation and neurogenesis; regulates transcription and chromatin remodeling activities required during neurogenesis (By similarity).[UniProtKB/Swiss-Prot Function]

shRNA Design:

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).