

Product datasheet for TL518003V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Flvcr1 Mouse shRNA Lentiviral Particle (Locus ID 226844)

Product data:

Product Type: shRNA Lentiviral Particles

Product Name: Flvcr1 Mouse shRNA Lentiviral Particle (Locus ID 226844)

Locus ID: 226844

Synonyms: 9630055N22Rik; FLVCR; Flvcr1; Mfsd7b

Vector: pGFP-C-shLenti (TR30023)

Format: Lentiviral particles

Components: Mfsd7b - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble

control), 0.5 ml each, >10^7 TU/ml.

RefSeq: <u>BC024751</u>, <u>NM 001081259</u>, <u>NM 001313747</u>, <u>NM 198025</u>, <u>NM 001081259.1</u>, <u>BC024751.1</u>,

BC157992, BC010797, BC156476, BC172689

UniProt ID: B2RXV4

Summary: Isoform 1: Heme transporter that exports cytoplasmic heme. It can also export

coproporphyrin and protoporphyrin IX, which are both intermediate products in the heme biosynthetic pathway. Does not export bilirubin. Heme export depends on the presence of HPX and is required to maintain intracellular free heme balance, protecting cells from heme toxicity. Heme export provides protection from heme or ferrous iron toxicities in liver, brain, sensory neurons and during erythtopoiesis, a process in which heme synthesis intensifies.

Causes susceptibility to FeLV-C in vitro.[UniProtKB/Swiss-Prot Function]

shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To

be certain that your variant of interest is targeted, please contact <u>techsupport@origene.com</u>. If you need a special design or shRNA sequence, please utilize our <u>custom shRNA service</u>.



Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).