

Product datasheet for **TL518019**

Clasp1 Mouse shRNA Plasmid (Locus ID 76707)

Product data:

| | |
|---------------------------|--|
| Product Type: | shRNA Plasmids |
| Product Name: | Clasp1 Mouse shRNA Plasmid (Locus ID 76707) |
| Locus ID: | 76707 |
| Synonyms: | 1700030C23Rik; 5730583A19Rik; B130045P17Rik; mKIAA0622 |
| Vector: | pGFP-C-shLenti (TR30023) |
| E. coli Selection: | Chloramphenicol (34 ug/ml) |
| Mammalian Cell Selection: | Puromycin |
| Format: | Lentiviral plasmids |
| Components: | Clasp1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 76707). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free. |
| RefSeq: | BC094432 , NM_001081276 , NM_001293300 , NM_001293301 , NM_029709 , NM_177548 , NM_001359328 , NM_001359329 , NM_001359330 , NM_001359331 , NM_001359334 , NM_001359335 , NM_001359336 , NM_001359337 , NM_001359338 , NM_001359339 , NM_001359340 , NM_001359341 , NM_001359342 , NM_001359343 , NM_001359345 , NM_001359346 , NM_001359347 , NM_029709.1 , NM_029709.2 , NM_177548.1 , NM_177548.2 , NM_001081276.1 , NM_001293301.1 , NM_001293300.1 , BC035533 , BC039216 , BC057312 , BC075708 , BC094554 , BC157938 , BC158064 |
| UniProt ID: | Q80TV8 |
| Summary: | Microtubule plus-end tracking protein that promotes the stabilization of dynamic microtubules. Involved in the nucleation of noncentrosomal microtubules originating from the trans-Golgi network (TGN). Required for the polarization of the cytoplasmic microtubule arrays in migrating cells towards the leading edge of the cell. May act at the cell cortex to enhance the frequency of rescue of depolymerizing microtubules by attaching their plus-ends to cortical platforms composed of ERC1 and PHLDB2. This cortical microtubule stabilizing activity is regulated at least in part by phosphatidylinositol 3-kinase signaling. Also performs a similar stabilizing function at the kinetochore which is essential for the bipolar alignment of chromosomes on the mitotic spindle (By similarity).[UniProtKB/Swiss-Prot Function] |



[View online »](#)

- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).