

Product datasheet for **TL518387**

Myof Mouse shRNA Plasmid (Locus ID 226101)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Myof Mouse shRNA Plasmid (Locus ID 226101)
Locus ID:	226101
Synonyms:	2310004N10Rik; 2310051D19Rik; E030042N20Rik; Fer1; Fer1l3
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Myof - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 226101). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001099634 , NM_001302140 , NM_001099634.1 , NM_001302140.1 , BC025449 , BC025649 , BC044825 , BC051217 , BC055953 , BC150790 , NM_177035
UniProt ID:	Q69ZN7
Summary:	The protein encoded by this gene is a member of the ferlin family of proteins, which have been implicated in fusion events in muscle tissue. Members of this family have a carboxy-terminal single pass transmembrane domain and multiple C2 domains, which bind negatively charged phospholipids in the presence of calcium ions. This gene is expressed at high levels in myoblasts and upregulated in damaged skeletal muscle. Mice deficient in this protein display defects in myoblast fusion, muscle regeneration, and angiogenesis. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Oct 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).