

Product datasheet for **TL518798**

Cecr2 Mouse shRNA Plasmid (Locus ID 330409)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Cecr2 Mouse shRNA Plasmid (Locus ID 330409)
Locus ID:	330409
Synonyms:	2610101O16Rik; 2810409N01Rik; Gtl4; mKIAA1740
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Cecr2 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 330409). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001128151 , NM_001128151.1 , BC060152 , NM_001368706 , NM_001128151.2
UniProt ID:	E9Q2Z1
Summary:	Chromatin reader component of histone-modifying complexes, such as the CERF (CECR2-containing-remodeling factor) complex and ISWI-type complex (By similarity). It thereby plays a role in various processes during development: required during embryogenesis for neural tube closure and inner ear development (PubMed:15640247, PubMed:20589882, PubMed:21246654). In adults, required for spermatogenesis, via the formation of ISWI-type chromatin complexes (PubMed:22154806). In histone-modifying complexes, CECR2 recognizes and binds acylated histones: binds histones that are acetylated and/or butyrylated (By similarity). May also be involved through its interaction with LRPPRC in the integration of cytoskeletal network with vesicular trafficking, nucleocytosolic shuttling, transcription, chromosome remodeling and cytokinesis (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).