

## Product datasheet for **TL519052**

### Flna Mouse shRNA Plasmid (Locus ID 192176)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Flna Mouse shRNA Plasmid (Locus ID 192176)
Locus ID:	192176
Synonyms:	ABP-280; Dilp2; F730004A14Rik; filamin-1; Fln1; GENA 379
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Flna - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 192176). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001290421</a> , <a href="#">NM_010227</a> , <a href="#">NM_010227.1</a> , <a href="#">NM_010227.2</a> , <a href="#">NM_010227.3</a> , <a href="#">NM_001290421.1</a> , <a href="#">BC151024</a> , <a href="#">BC004061</a> , <a href="#">BC038478</a> , <a href="#">BC039208</a> , <a href="#">BC054432</a> , <a href="#">BC099506</a>
Summary:	Actin binding protein that promotes orthogonal branching of actin filaments and links actin filaments to membrane glycoproteins. Anchors various transmembrane proteins to the actin cytoskeleton and serves as a scaffold for a wide range of cytoplasmic signaling proteins (By similarity). Interaction with FLNA may allow neuroblast migration from the ventricular zone into the cortical plate. Tethers cell surface-localized furin, modulates its rate of internalization and directs its intracellular trafficking. Involved in ciliogenesis. Plays a role in cell-cell contacts and adherens junctions during the development of blood vessels, heart and brain organs (PubMed:17172441). Plays a role in platelets morphology through interaction with SYK that regulates ITAM- and ITAM-like-containing receptor signaling, resulting in by platelet cytoskeleton organization maintenance (PubMed:20713593).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).