

Product datasheet for **TL519213**

Qk Mouse shRNA Plasmid (Locus ID 19317)

Product data:

Product Type:	shRNA Plasmids
Locus ID:	19317
Synonyms:	1110003F05Rik; 1500005P18; l(17)-1Wis; l17Wisl; Qkl
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Qk - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector (Gene ID = 19317). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC053426 , BC056346 , NM_001159516 , NM_001159517 , NM_021881 , NM_001159517.1 , NM_021881.1 , NM_021881.2 , NM_001159516.1 , BC053426.1
UniProt ID:	Q9QYS9
Summary:	RNA-binding protein that plays a central role in myelinization (PubMed:10864952, PubMed:11917126). Also required for visceral endoderm function and blood vessel development (PubMed:11892011, PubMed:16470614). Binds to the 5'-NACUAAAY-N(1,20)-UAAY-3' RNA core sequence (PubMed:16041388). Acts by regulating pre-mRNA splicing, mRNA export, mRNA stability and protein translation, as well as cellular processes including apoptosis, cell cycle, glial cell fate and development (PubMed:10535969, PubMed:12467586, PubMed:11297509, PubMed:11917126, PubMed:15568022). Required to protect and promote stability of mRNAs such as MBP and CDKN1B which promotes oligodendrocyte differentiation (PubMed:10535969, PubMed:15568022). Participates in mRNA transport by regulating the nuclear export of MBP mRNA (PubMed:12467586). May also play a role in smooth muscle development (PubMed:14706070).[UniProtKB/Swiss-Prot Function]



shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

Performance Guaranteed: OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).