

Product datasheet for **TL519327V**

Ufm1 Mouse shRNA Lentiviral Particle (Locus ID 67890)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Ufm1 Mouse shRNA Lentiviral Particle (Locus ID 67890)
Locus ID:	67890
Synonyms:	1810045K17Rik; AI132708; AI463323; ENSMUSG00000074598; Gm10726
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Ufm1 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	BC061065 , NM_026435 , NM_026435.1 , NM_026435.2 , NM_026435.3 , NM_026435.4 , NM_026435.5 , BC033280 , BC051554
UniProt ID:	P61961
Summary:	Ubiquitin-like modifier which can be covalently attached via an isopeptide bond to substrate proteins as a monomer or a lysine-linked polymer (PubMed:21494687). The so-called ufmylation, requires the UFM1-activating E1 enzyme UBA5, the UFM1-conjugating E2 enzyme UFC1, and the UFM1-ligase E3 enzyme UFL1. This post-translational modification on lysine residues of proteins may play a crucial role in a number of cellular processes. TRIP4 ufmylation may for instance play a role in nuclear receptors-mediated transcription (By similarity). Other substrates may include DDRGK1 with which it may play a role in the cellular response to endoplasmic reticulum stress.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).