

Product datasheet for **TL519417**

Cpvl Mouse shRNA Plasmid (Locus ID 71287)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Cpvl Mouse shRNA Plasmid (Locus ID 71287)
Locus ID:	71287
Synonyms:	4933436L16Rik; HVLP
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Cpvl - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 71287). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001289713 , NM_001289714 , NM_027749 , NM_027749.1 , NM_027749.2 , NM_001289714.1 , NM_001289713.1 , BC137839 , BC144966
UniProt ID:	Q9D3S9
Summary:	This gene encodes a member of the serine carboxypeptidase family of proteases that cleave amino acids from the C-terminus of a protein substrate. The human ortholog of this gene, where it was first characterized, was found to be upregulated during the maturation of monocytes to macrophages. The encoded protein may be involved in antigen processing, digestion of phagocytosed proteins in the lysosome and lamellipodium formation. Disruption of this gene in mice was found to cause embryonic lethality. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Aug 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).