

Product datasheet for **TL519430**

Zswim8 Mouse shRNA Plasmid (Locus ID 268721)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Zswim8 Mouse shRNA Plasmid (Locus ID 268721)
Locus ID:	268721
Synonyms:	2310021P13Rik; 4832404P21Rik; Kiaa0913; mKIAA0913
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Zswim8 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 268721). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001252081 , NM_001252082 , NM_027996 , NM_175544 , NM_027996.1 , NM_027996.2 , NM_027996.3 , NM_001252082.1 , NM_001252081.1 , BC151173 , BC023935 , BC026504 , BC027155 , BC030914 , BC037627 , BC049362 , BC058666 , BC059058 , BC085161 , BC151046 , BC151056 , BC172103
UniProt ID:	Q3UHH1
Summary:	Substrate recognition component of a SCF-like E3 ubiquitin-protein ligase complex that promotes target-directed microRNA degradation (TDMD), a process that mediates degradation of microRNAs (miRNAs) (PubMed:33184237). The SCF-like E3 ubiquitin-protein ligase complex acts by catalyzing ubiquitination and subsequent degradation of AGO proteins (AGO1, AGO2, AGO3 and/or AGO4), thereby exposing miRNAs for degradation (By similarity). Specifically recognizes and binds AGO proteins when they are engaged with a TDMD target (By similarity). May also acts as a regulator of axon guidance: specifically recognizes misfolded ROBO3 and promotes its ubiquitination and subsequent degradation (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).